

Afiscooter S

The ultimate outdoor Heavy-duty mobility scooter

The Afiscooter S with its spacious design maximize the protection and safety of the driver, with its brilliant ergonomics and with its excellent comfort of driving is one of the leading mobility scooters in its category.

This luxurious vehicle that has a proven record of high quality and long-lasting reliability can travel comfortably and safely in various types of terrains as rough off-road terrain (using golf tires) and extreme weather conditions (when using the unique hardtop canopy).

With single or wide seat, 3 wheel model (maximum maneuverability) or 4 wheel model (the ultimate in road presence and stability), either one you choose would fit you perfectly.

The Afiscooter SOLAR the first solar mobility scooter using sustainable energy, allowing up to 30% charging. The Afiscooter SOLAR allowing day use on solar power.

Experience safer driving with the CSA - Curve Speed Adaptation system.

Our innovative system automatically reduces speed on curves & sharp turns, ensuring your safety on the road.



Pleasurable driving experience



Designed for safety



Ease of operation



Top quality & Reliability



Highlights

- New rubber (pu) shock absorbing front bumper to maximize the protection and safety of the driver.
- An ergonomic tiller with a bright, easy to read and improved LCD display.
- A luxurious orthopedic rotatable seat with a high backrest that provides comfort while driving and ease boarding (offered in three widths: 20", 22" & 24").
- Powerful, long-lasting, energy efficient led headlights.
- The option of wide gt tires, for extra traction and stability together with the high chassis clearance give the afiscooter s off road performance.
- Optional hardtop canopy with ergonomic design to ensure comfortable ride in all weather conditions.
- Extremely spacious.



Afikim electric vehicles is a leader in developing and manufacturing premium, electrically-powered mobility scooters for more than 40 years. The well-known Afiscooter product line is designed with driver safety in mind first and foremost, alongside attention to ergonomics, comfort and reliability of the vehicle. Afikim electric vehicles has established a worldwide

distribution network that offers local technical support and follows the highest standards of service. Our company meets highest standards of quality and European ISO 9001, CE and EN 12184 certifications, as well as US safety standards with FDA approval.

Offering a safe and pleasurable driving experience on and off the road



Features



Specification

| | | Afiscooter S3 | | Afiscooter S4 | |
|---------------------|------------------|--|---------------|--------------------|---------------|
| | | Single | Wide Seat | Single | Wide Seat |
| Overall Length | | 155cm (61") | | 165cm (65") | |
| Overall Width | Standard Wheels | 68cm (27") | 96cm (37") | 68cm (27") | 96cm (37") |
| | Golf Tires | 88cm (35") | 96cm (37") | 88cm (35") | 96cm (37") |
| Overall Height | (-) Canopy | 125cm (49") | | 125cm (49") | |
| | (+) Canopy | 164cm (65") | | 164cm (65") | |
| Ground Clearance | | 12.5cm (5") | | 12.5cm (5") | |
| Weight | (+) Batteries | 153kg (337lb) | 163kg (359lb) | 155kg (342lb) | 165kg (364lb) |
| | (-) Batteries | 106kg (234lb) | 116kg (256lb) | 108kg (238lb) | 118kg (260lb) |
| | (+) Canopy | +15Kg (33lb) | | | |
| | (+) Wide Wheels | +12kg (+26lb) | | | |
| Weight Capacity | Incl. Driver | S: 200kg (500lb) S-HD: 250kg (600lb) | | | |
| Travel Range | (+) Batteries | Up to 60km (37m) | | | |
| Maximum Speed | | 10-18Kph (6.2-11mph)* | | | |
| Max Climbing Angle | | 8° (15%) 15Kph, 11° (20%) 12Kph | | | |
| Turning Radius | | 143cm (56") | | 228cm (90") | |
| Seating Width | | 51cm (20") | 85cm (33") | 51cm (20") | 85cm (33") |
| Tire Type | | Pneumatic (1.7bar) | | Pneumatic (1.7bar) | |
| Tire Size | Front Tyres | 3-10" | | 4.10/3.50-6" | |
| | Rear Tyres | 3.50-10" | | 3.50-10" | |
| Battery Type | Two | 12V 105Ah, Max. | | | |
| Power Unit | Rear Wheel Drive | 24V DC motor 1,400W | | | |
| Controller Type | | SD 140A | | | |
| Charger Type | | 8Amp | | | |
| Options/Accessories | | Left speed lever. Pedal operation Canopy Winter light sides close (Rain Cover) Wide seats 22" & 24" Golf Bag holder Crutch Holder Golf Tires | | | |

* Speed limit varies between the different countries according to the local regulation of each county.

** U.S model | Tested by manufacturer, based on actual parameters | The above specifications subject to change without notice