

# Afiscooter

The ultimate outdoor  
Heavy-duty mobility scooter



The **Afiscooter S** with its spacious design maximize the protection and safety of the driver, with its brilliant ergonomics and with its excellent comfort of driving is **one of the leading mobility scooters in its category.**

This luxurious vehicle that has a proven record of high quality and long-lasting reliability can travel comfortably and safely in various types of terrains as rough off-road terrain (using golf tires) and extreme weather conditions (when using the unique hardtop canopy).

With single or wide seat, 3 wheel model (maximum maneuverability) or 4 wheel model (the ultimate in road presence and stability), either one you choose would fit you perfectly.

## Highlights

- New rubber (pu) shock absorbing front bumper to maximize the protection and safety of the driver.
- An ergonomic tiller with a bright, easy to read and improved LCD display.
- A luxurious orthopedic rotatable seat with a high backrest that provides comfort while driving and ease boarding (offered in three widths: 20", 22" & 24").
- Powerful, long-lasting, energy efficient led headlights.
- The option of wide gt tires, for extra traction and stability together with the high chassis clearance give the afiscooter s off road performance.
- Optional hardtop canopy with ergonomic design to ensure comfortable ride in all weather conditions.
- Extremely spacious.



Pleasurable  
driving experience



Design for  
safety



Ease of  
operation



Top quality  
& Reliability



**Afikim electric vehicles** is a leader in developing and manufacturing premium, electrically-powered mobility scooters for more than 40 years. The well-known Afiscooter product line is designed with driver safety in mind first and foremost, alongside attention to ergonomics, comfort and reliability of the vehicle. Afikim electric vehicles has established a worldwide distribution network that offers local technical support and follows the highest standards of service.

Our company meets highest standards of quality and European ISO 9001, CE and EN 12184 certifications, as well as US safety standards with FDA approval.

Offering a safe and pleasurable driving experience on and off the road



## Features



## Specification

		Afiscooter S3		Afiscooter S4	
		Single	Wide Seat	Single	Wide Seat
Overall Length		155cm (61")		165cm (65")	
Overall Width	Standard Wheels	68cm (27")	96cm (37")	68cm (27")	96cm (37")
	Golf Tires	88cm (35")	96cm (37")	88cm (35")	96cm (37")
Overall Height	(-) Canopy	125cm (49")		125cm (49")	
	(+) Canopy	164cm (65")		164cm (65")	
Ground Clearance		12.5cm (5")		12.5cm (5")	
Weight	(+) Batteries	153kg (337lb)	163kg (359lb)	155kg (342lb)	165kg (364lb)
	(-) Batteries	106kg (234lb)	116kg (256lb)	108kg (238lb)	118kg (260lb)
	(+) Canopy	+15kg (33lb)			
	(+) Wide Wheels	+12kg (+26lb)			
Weight Capacity	Incl. Driver	S: 200kg (450lb)   S-HD: 250kg (550lb)			
Travel Range	(+) Batteries	Up to 60km (37m)			
Maximum Speed		10-18Kph (6.2-11mph)*			
Max Climbing Angle		8° (15%) 15Kph, 11° (20%) 12Kph			
Turning Radius		143cm (56")		228cm (90")	
Seating Width		51cm (20")	85cm (33")	51cm (20")	85cm (33")
Tire Type		Pneumatic (1.7bar)		Pneumatic (1.7bar)	
Tire Size	Front Tyres	3-10"		4.10/3.50-6"	
	Rear Tyres	3.50-10"		3.50-10"	
Battery Type	Two	12V 105Ah, Max.			
Power Unit	Rear Wheel Drive	24V DC motor   1,400W			
Controller Type		SD 140A			
Charger Type		8Amp			
Options/Accessories		Left speed lever.   Pedal operation   Canopy   Winter light sides close (Rain Cover)   Wide seats 22" & 24"   Golf Bag holder   Crutch Holder   Golf Tires			

\* Speed limit varies between the different countries according to the local regulation of each country.

\*\* U.S. model | Tested by manufacturer, based on actual parameters | The above specifications subject to change without notice