

The Afiscooter

An icon in the world of personal electric mobility - experience you can count on



A balance of simplicity, reliability and classic design marks the **Afiscooter SE** as a world leading personal electric vehicle for over 40 years.

The strong steel chassis and rear power unit combine strength with a long driving range, designed for stable and safe ride on any terrain.



Highlights

- Reliable electronic command system.
- Comfortable seat with adjustable armrest.
- Fully electric- provides a silent, environmental friendly and highly efficient maintenance vehicle.
- Three wheel model for maximum maneuverability.
- Solid large rear basket for storage.



Design for safety



Ease of operation



Icon in the world of personal mobility



Afikim Electric Vehicles is a leader in developing and manufacturing premium, electrically-powered mobility scooters for more than 40 years. The well-known Afiscooter product line is designed with driver safety in mind first and foremost, alongside attention to ergonomics, comfort and reliability of the vehicle. Afikim Electric Vehicles has established a worldwide distribution network that offers local technical support and follows the highest standards of services.

Our company meets highest standards of quality and European ISO 9001, CE and EN 12184 certifications, as well as US safety standards with FDA approval.



Features



Specification

		Afiscooter SE	
		Single	Dual Seat
Overall Length		158cm (62")	166 cm (65")
Overall Width	Standard Wheels	73cm (29")	91 cm (36")
Overall Height	(-) Canopy	124cm (49")	101 cm (40")
	(+) Canopy	202cm (79")	180cm (70")
Ground Clearance		12.5cm (5")	
Weight	(+) Batteries	146kg (322lb)	163kg (359lb)
	(-) batteries	99kg (218lb)	116kg (256lb)
Weight Capacity	Incl. Driver	200kg (500lb)	
Travel Range		45km (28m)	
Maximum speed		15km/h (9.3mph)	12km/h (7.5mph)
Max Climbing Angle		20%	
Turning Radius		115cm (45")	
Seating Width		46cm (18")	85cm (33")
Tire type		Pneumatic	
Tire size	Front Tyres	2.75x10"	
	Rear Tyres	2.75x12"	
Battery type	two	12V 73Ah	
Power unit	Rear wheel drive	24V DC motor 1,300W	
Charger type		8Amp	
Color		Black & Chrome	

* Speed limit varies between the different countries according to the local regulation of each county.

** U.S model | Tested by manufacturer, based on actual parameters | The above specifications subject to change without notice