




Large brake bar



Easily foldable frame

Indoor Walker

KA382WHI

Features

- Lightweight aluminium frame
- Height adjustable handle for the perfect fit
- Ergonomically designed handle for comfort
- Large brake bar with plastic cover for improved grip and feel
- Easily foldable compact frame
- Basket for storage
- Equipped with a convenient tray for effortlessly carrying small items, food and drinks
- White frame with a black coloured tray as standard
- Wood patterned trays are available as an option

Specifications

Product Net Weight	5kg
Overall Height	880-920mm
Overall Width	600mm
Overall Length	670mm
Handle Width	450mm
Tray Size	W385 x D285mm
Tray Height (to floor)	720mm
Folded Dimension	L210 x W600 x H990mm
Castors	150mm
Max. User Weight	100kg
Tray Maximum Load	5kg
Basket Maximum Load	5kg
Main Frame Material	Aluminium

Please contact customerservice@kcare.com.au for more information.
www.kcare.com.au | 1300 783 783