

FOLD & GO

COMPACT FOLDING POWERCHAIR
P113













FEATURES

- Lightweight, foldable power wheel chair which folds easily in seconds for transportation or storage
- Folds in one easy step
- Compact size means great manoeuvrability in tight spaces
- Removable battery box with battery charge indicator & extra charge port for charging off board
- Adjustable height footplate
- Adjustable height backrest
- Intelligent, regenerative, electromagnetic braking
- Vibrant smart colours available, red, blue or grey

SPECIFICATION

Max. Range	11km / 7mi. 15km / 9.6mi. (optional)
Max. Speed	7.5kph / 4.7kph
Battery	1pc Lithium-ion battery 24V/11.4AH / 17.1AH (optional)
Motor	DC 24V, 70W (Rating)/210W (Max)
Gradient	3°
Controller	PG VR2 50A
Charger	2A
Caster	7" PU tire
Drive Wheel	8"(200x50) PU tire
Brake	electromagnetic brakes

DIMENSIONS (cm/in)

 Seat Width	 Seat Depth	 Seat Height from Ground	
46/18	42/16.5	52/20.5	
 Armrest Height	 Backrest Height	 Footrest Extension	
22/9	33/13	42/16.5 48/19	
 Overall Length	 Overall Width	 Overall Height	 Folded Dimension
97cm/38"	59/23	95-100/ 37.5-39	59x37x81/ 23x15x32 LxWxH

WEIGHT (kg/lb)

tll w/batt.	tll w/o batt.	capacity
23/50.5 23.9/52.5 (optional)	20.9/46	114/250

COMPACT FOLDING POWER WHEELCHAIR



ISO 13485
ISO 14001
rvv2018-1