

Maverick 14




Maverick 14

Features

- Mid-wheel drive design offers outstanding manoeuvrability in tight, compact spaces and the superior adjustable suspension design delivers a smooth and stable ride both indoors and outdoors, even over tough terrain
- Patented obstacle climbing suspension that works
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Swing away joy stick standard
- Height and width adjustable armrests on the new contoured seat. Height adjustable four post seat with reclining back rest
- 204kg weight capacity
- Available in a number of vibrant new colours – Red, Silver, Pearl Gold, Charcoal, NV Green, Kite Blue
- Other options available are electric seat lift, rear attendant controls, optional seating & leg rests



Dimensions cm / in



Overall Dimension : A x B x C	101 x 66 x 127cm / 40 x 26 x 50in
Seat Dimension : D x E x F	51 x 46 x 51cm / 20 x 18 x 20in
Seat Height from Ground : G	56-58cm / 22-23in
Seat Height from Deck : H	46-51 cm / 18-20in

Weight kg / lb

base	battery/ea	seat	total with battery	capacity
60/132	17/37	21/46	115kg / 253lb	205kg / 450lb

Technical Specs.

Max Speed	7.2kph / 4.5mph
Range	40km/ 25mi
Turn Radius	53cm/ 21"
Ground clearance	8cm/ 3.1"
Motor	DC 24V, 200W(Rating) / 575W(Max)
Controller	Dynamic Shark 50
Battery	12V/55-62AH x 2pcs
Charger	5A off-board
Gradient	8°
Caster	front: 8" x 2" PU tire ; rear: 8" PU tire
Drive Wheel	14" tire
Brake	Intelligent, regenerative, electromagnetic brakes