



Eco 745


Eco 745

Features

- Stylish eye catching design with attractive modern colour range
- Loop handlebars offer comfortable finger or thumb throttle operation
- Adjustable high back contoured seat with simple operation of slide and swivel mechanisms to allow maximum user comfort
- Infinite adjustment of steering tiller allowing everyone to easily find a comfortable position
- High mount charger socket means easy access when charging
- Removable front basket with carry handle
- Adjustable power saving LED headlight, indicators front and rear, rear tail lights and reverse warning
- Half speed switch fitted.



Dimensions cm/in



Overall Dimension A x B x C	120 x 61 x 99cm/47 x 24 x 39"
Seat Dimension : D x E x F	46 x 40 x 46cm/18 x 16 x 18"
Seat Height from Ground : G	53-58cm/21-23"
Seat Height from Deck : H	40-47cm/16-18.5"
Wheel Base : I	85cm/33.5"

Weight kg/lb

base	battery/ea	seat	total with battery	capacity
51/112	10/22	15/33	86/190	136/300

Technical Specs

Max Speed	10kph/6.25mph
Range	20-30km/13-20mi
Turn Radius	157cm/62"
Ground clearance	9.5cm/3.74"
Motor	DC 24V, 210W (Rating)/715W (Max)
Controller	S-Drive 70A
Battery	12V/34A/H x 2pcs
Charger	5A off-board
Gradient	6°
Front Wheel	10" Pneumatic tire
Rear Wheel	10" Pneumatic tire
Brake	intelligent, regenerative, electromagnetic brakes