

Boston



Features:

- ✓ Designed with front & rear suspension to give a more comfortable ride even on flat surfaces
- ✓ Designed to pull apart with ease for transportation or storage
- ✓ Modular design pulls apart easily & can fit into the car boot – tiller folds down and seat lifts off
- ✓ Comfortable seat is height adjustable to help you find the right driving position and swivels 360 degrees for easy mounting and dismounting
- ✓ Solid, puncture proof tyres are easy to keep clean and maintain
- ✓ Two 12V 12ah batteries provide power to travel around the local shopping centre and other shorter journeys
- ✓ Horn is useful for alerting nearby pedestrians to your presence when driving in crowded areas
- ✓ Anti-tippers help to minimize risk of tipping on corners or bumpy ground
- ✓ Available in a choice of colours – red, sporty graphite grey and bright white

Tech Specs

Dimensions	L102cm x W61cm x H90cm
Scooter Total Weight	48kg
Battery Pack Weight	10kg
Rear Assy Weight	10.5kg
Max Speed	8kph
Max Range	10-15km
Max Weight Capacity	136kg
Turning Circle	120cm
Ground Clearance	8cm
Motor	220w
Batteries	2x12V 12AH
Charger	24V / 2A



ROTHCARE

