

Lite Ryder - Lightweight Folding Powerchair



Features:

- ✔ Super Lightweight carbon fibre frame
- ✔ Folds in one easy step for transportation
- ✔ Comfortable padded seat & backrest
- ✔ Lightweight removable Lithium battery
- ✔ Compact easy to use joystick
- ✔ Complete chair weighs less than 13kg with the battery
- ✔ USB fitting in Joystick & in removable battery pack
- ✔ Battery pack is fitted with a charging port for offboard charging

Australian Standard compliant AS/NZS 3695.2:2019

Tech Specs



Overall width	53cm
Overall length	90cm
Overall height	86cm
Seat width	43cm
Seat depth	40cm
Seat height from ground	46cm
Backrest height	40cm
Speed	6kph - 10kph
Turning radius	90cm
Gradient	3 degrees
Range	10km
Battery	24V / 10Ah Lithium-ion
Motor	DC 24V, 150W brushless DC
Controller	Shanghai Micon
Brake	Electromagnetic
Front wheel	6"
Rear wheel	10.6"
Weight capacity	100kg
Weight w/out battery	11.85kg
Weight with battery	13kg

