



KCare

Electric Tilt
in Space


Mobile Shower Commode

**COMFORT
SUPPORT AND
SAFETY**

KCare
Healthcare Solutions



IPX6 Linak powered electric actuator tilt mechanism



3x 125 mm twin locking castors with 1x 125mm directional locking castor



Adjustable height swing-away foot rests with padded calf rests

Electric Tilt in Space Mobile Shower Commode

KA11506

KA115065

Features

- IPX6 Linak powered electric actuator tilt mechanism
- Control of tilting action is performed by attendant operated control unit
- Extra-long handset cable allows the carer to perform tilting action whilst engaging with the user
- Tilt from 0 to 25 degrees from vertical
- Curved padded backrest
- Removable arms with padded arm rests
- Adjustable height swing-away foot rests with padded calf rests
- Adjustable padded head rest
- Integrated push handle
- 3x 125 mm twin locking castors with 1x 125mm directional locking castor make the unit easy to manoeuvre and secure when transferring
- Optional open or closed front seat in standard or pressure relieving densities
- Maximum user weight of 300kg

Specifications

KA11506

Seat width	600mm
Seat depth	470mm
Floor to top of seat frame	550mm
Backrest height	630mm
Tilt angle	0-25°
Overall width	710mm
Overall depth	1140mm
Overall height	1250mm
Max user weight	300kg

KA115065

Seat width	650mm
Seat depth	480mm
Floor to top of seat frame	550mm
Backrest height	630mm
Tilt angle	0-25°
Overall width	760mm
Overall depth	1140mm
Overall height	1250mm
Max user weight	300kg

Please contact customerservice@kcare.com.au or call **1300 783 783** to find out more information on specific models in each range.

836a Fifteenth Street Mildura VIC 3500