

General Advice: These instructions are to be used in conjunction with the information your Health Professional has discussed with you. All equipment is to be used following the manufacturers weight recommendations. If you are unsure on how to use or fit a piece of equipment, please contact your Allied Health Professional.

Day Chair Use

- A day chair is a type of chair that is designed for use during the day, typically in a living room or common area of a home
- Day chairs can be used for reading, watching television, or other leisure activities. They are available in different materials, sizes and are height adjustable to suit the functionality of the user

Installation

1. Your day chair should come assembled, but depending on height requirements you may need to adjust the height of your chair
2. Following the user manual instructions, raise or lower your chair to suit the user. Your Allied Health Professional or a member of the Country Care Group service team can assist in ensuring the height is suitable

How to Use a Day Chair

1. Adjust the legs on the chair using the steel buttons to slide the internal leg piece up or down to the appropriate height. The users feet should be in contact with the floor
 2. Some day chairs have an adjustable backrest to make the seat depth longer or shorter, this can be adjusted with screws underneath the chair
 3. Remember to clean and maintain the chair as per the manufacturer's instructions. This can help to ensure that the chair lasts for a long time
- **Note:** Day chairs may vary in design, comfort and functionality, so you should familiarise yourself with the unique features of your day chair and adjust accordingly

Precautions & Safety

- Ensure that the day chair is stable on the floor & cannot slide
- Ensure there is adequate space around the day chair to enable safe movement on, off and around the chair

Maintenance & Cleaning

- Please check the product label / user manual for specific cleaning instructions
- If the equipment appears damaged, **DO NOT USE**. It may require replacing



Complementary Assistive Devices



Sensor Light



Cushion



Reachers



Foot Stool



Over Chair Table



Kitchen Tray Mobile



Chair Protectors



One Touch Can Opener



Tipping Kettle