

General Advice: These instructions are to be used in conjunction with the information your Health Professional has discussed with you. All equipment is to be used following the manufacturers weight recommendations. If you are unsure on how to use or fit a piece of equipment, please contact your Allied Health Professional.

Swivel Bath Seat Use

- A swivel bath seat is designed to assist people with limited mobility or disability in bathing. This device is a bath seat that rotates/swivels, allowing the person to be transferred in and out of the bathtub with ease and safety

Installation

- Installing a swivel bather requires proper training and instruction, so it is important to consult with a healthcare professional and read the manufacturer's instructions and user guide before installing the device

Here are some general steps on how to install a swivel bath seat:

1. Measure the inside of your bathtub to ensure that the swivel bather will fit properly
2. Place the swivel bath seat base (without the seat) into the bathtub, making sure that it is centered and level
3. Secure the base to the bathtub using the positioning screws located on both sides. The black rubber ends should be screwed against the side of the bath, which should prevent the base from moving
4. Place the seat on to the socket of the base. Check the swivel bath seat is securely attached and that it swivels and rotates smoothly by pulling up the lever
5. If you require increased height to the swivel bath seat an 100mm raiser can be installed (ask supplier if this is required)
6. Make sure that all the components are securely fastened and tightened

How to Use a Swivel Bath Seat

1. Before using the swivel bath seat, make sure that it is securely attached and adjusted to a comfortable and safe position for the person who will be using it
2. Sit on the swivel bath seat
3. Use the lever on the side to rotate into the bathtub, legs will need to be lifted over and into the bath during rotation either independently by the user or with assistance
4. Once the user is finished bathing, use the swivel function to rotate to the original position, lifting legs back out and over the bathtub
5. It is not recommended at any time to stand from the swivel bath seat while in the bath

Precautions & Safety

- Ensure that the seat is locked into position before sitting or standing
- Ensure that items required for showering are in easy reach when seated

Maintenance & Cleaning

- Please check the product label /user manual for specific cleaning instructions
- If the equipment appears damaged, **DO NOT USE.** It may require replacing

